


## Powerful System for Quantitative Analysis of Thin-Layer Chromatograms



Application

Development

Derivatization

Documentation

For quantitative determination of samples, the **HPTLC-Densitometer CD60** converts the spots/bands of the single substances into a chromatogram curve. It measures the absorbance or fluorescence of separated compounds in transmission or reflection mode. The **HPTLC-Densitometer CD60** is controlled by ProQuant software which also enables quantitative evaluation of the generated data.

### Key Aspects of HPTLC-Densitometer CD60

- Absorbance or fluorescence measurement
- Remission and transmission mode
- For objects up to 265 x 200 x 4 mm
- Automatic starting and switching of filters and lamps
- Rapid data collection and evaluation
- Recording spectra from 190 - 900 nm
- Software-controlled by ProQuant
- Ease of operation
- Reproducible and reliable results
- Meeting the requirements of GMP/GLP



The **HPTLC-Densitometer CD60** works within a spectral range of 190 - 900 nm. This is provided by three light sources: a deuterium lamp (190 to 340 nm), a halogen lamp (340 to 900 nm) as well as a mercury lamp. Once the wavelength is selected, the densitometer will automatically start to scan the entire plate. It measures the absorbance

or fluorescence reflected or transmitted by each sample. This will be stored in the software in the form of peak tables. These tables consist of Rf values and area of each spot. Therefore, you can carry forward the quantitative evaluation of the generated densitometric data by ProQuant software.

## Measurement and Evaluation

### Method Types



- Method for chromatogram
- Method for multi-wavelength scan
- Method for spectrum

### Recording Modes

- Remission and transmission
- Absorbance or fluorescence
- Linear and Meander scan
- Two-wavelength measurement
- Multi-wavelength measurement

### Results

- Peak lists
- Results for sample and standard
- Automatic integration with manual correction facility
- Linear, polynomial or Michaelis-Menten function




## Technical Parameters of HPTLC-Densitometer CD60

|                        |                                |
|------------------------|--------------------------------|
| Object size            | Up to 265 x 200 x 4 mm         |
| Spectral range         | 190 - 900 nm                   |
| Filters                | 370, 420, 450, 550, orange, UV |
| Max. scan length       | 5 to 195 mm                    |
| Max. scan width        | 5 to 260 mm                    |
| Slit width             | 0.4 to 10 mm                   |
| Slit height            | 0.02 to 2 mm                   |
| Dimensions (W x H x D) | 730 x 550 x 300 mm             |
| Weight                 | 30 kg                          |

## Ordering Information

|           |  |
|-----------|--|
| BS131.800 | <b>HPTLC-Densitometer CD60</b> , 230 V, incl. interface box, software ProQuant |
| BS131.801 | <b>HPTLC-Densitometer CD60</b> , 110 V, incl. interface box, software ProQuant |
| BS131.816 | Software Provalid, program for automatic validation                            |
| BS131.830 | Software Spectra Calc, program for compilation of spectra libraries            |
| BS131.825 | IQ/OQ documents for HPTLC-Densitometer CD60                                    |

## Distributor



biostep GmbH  
 Innere Gewerbestraße 7  
 09235 Burkhardttsdorf, Germany  
 T +49 3721 39050, info@biostep.de