




Atomic Absorption Spectrophotometer

Fully Automatic Combined Flame and Graphite Furnace System

Model AAS-3900



Atomic Absorption Spectrophotometer AAS-3900 is a superior instrument for the research laboratory, and is an advanced and affordable system with high sensitivity that generates accurate and reproducible measurements. The AAS-3900 spectrophotometer is accurate, reliable, and is an exceptional value. With its built-in, computer-controlled Air/Acetylene flame, titanium alloy burner and high-efficiency glass nebulizer design, the system provides optimal and reproducible results for micro and macro samples with high resolution.

Atomic Absorption Spectrophotometer AAS-3900 has a **powerful built-in software** which permits this instrument to be linked to a computer and a printer to display the photometric and spectral data on the PC monitor. **Atomic Absorption Spectrophotometer AAS-3900's** enhanced transmission and full reflection makes this atomic absorption spectrophotometer highly effective and reduces noise.

One of its advantages is its accurate wavelength, ease of operation, versatile software applications, and effortless optional accessory installation. This instrument provides for automatic switching between the flame and graphite furnace, and is widely used for analyzing samples for **Agricultural, Food, Geological, Clinical, Metal, Petrochemical, Environmental, Mining and Pharmaceutical applications.**

It is easy to manipulate, and is fully automated, allowing for automatic adjustment of the lamp current and position, the burner head position, the negative voltage, and the gas flow. Safety is our primary concern, and the **Atomic Absorption Spectrophotometer AAS-3900** allows for constant monitoring of the burner head, the flame, the ignition, air pressure, and drain status, to ensure the optimum functioning of the instrument.

Atomic Absorption Spectrophotometer AAS-3900 has a highly effective nebulizer, the sensitivity of the Cu $5\mu/\text{ml}$ > 0.850 Abs. Labomed, Inc. is certified by ISO-9001-2013, has CE Conformity and is FDA Licensed.

Features

FEATURES AND FUNCTIONS:

The instrument has a motorized 8 hollow cathode lamp turret which allows the automatic positioning and optimization of each hollow cathode lamp by the software. The control of the gas flows for the fuel gas (C_2H_2) of the burner is also carried out directly from the software, thus allowing optimization of the instrument for the best analytical parameters for a selected analysis.

Two methods of background correction are available. The first utilizes a Deuterium Arc lamp and the second is the proven method of Self Reversal.

High precision minimal optics ensures maximum light throughput to the computer controlled Czerny-Turner Monochromator.

The location of the wavelength and peak selection is automatically controlled from the software.

The spectral bandwidth is automated and is available with a choice of five slit sizes.

The electronic parameters for the photomultiplier tube detector, the hollow cathode lamp current and the balancing of the absorbance and background energies are controlled from the software.

The ignition of the flame is computer controlled and the various safety interlocks offer a very safe operating system.

SAFETY:

The flame conditions are continuously monitored and should the flow rates change, an audible alarm sounds.

The pressure of the support gas (oxidant) is monitored constantly. If the pressure changes then the flow of the fuel gas will be stopped and the flame will be safely extinguished.

A sensor monitors the level of liquid in the drain and will prevent ignition if too slow. The flame will also be extinguished if the level of liquid in the drain changes significantly.

A flame sensor monitors the flame and safely turns off the gas flow to the burner if the flame suddenly extinguishes.

The burner is identified by a switch making it impossible to light without the burner being fitted.

An emergency flame off button is installed in case a problem is observed. The flame can be extinguished safely.



Atomic Absorption Spectrophotometer

Fully Automatic Combined Flame and Graphite Furnace System

Model AAS-3900

INCLUDED ACCESSORIES



Motorized 8 hollow cathode lamp turret accessory



Atomizer accessory



Graphite Furnace (interior)



Graphite Tube with Platform (close up)

OPTIONAL ACCESSORIES



HYDRIDE GENERATOR

A hydride generator is available for the determination of elements such as Arsenic, Selenium, Antimony, Tellurium and mercury at ultra low levels. The hydride generator is supplied with an absorption cell, and electrical absorption cell heater and controller and all necessary burner fittings.



GRAPHITE FURNACE AUTOSAMPLER

The graphite furnace auto-sampler system allows automatic update of calibration data, matrix modifications and automatic dilution of samples.



Atomic Absorption Spectrophotometer

Fully Automatic Combined Flame and Graphite Furnace System

Model AAS-3900

Optic System - Czerny-Turner (C-T) type Monochromator	
Wavelength Range	190nm - 900nm
Spectral Bandwidth	0.1nm, 0.2nm, 0.5nm, 1.0nm and 2.0nm (5 steps with automatic changeover).
Wavelength accuracy	+/- 0.3nm (automatic calculation)
Wavelength reproducibility	0.2nm
Resolution	Better than 0.3nm
Baseline stability	0.005A/30 min
Hollow cathode Lamp system	8 lamp simultaneously lit (one with pre-heat). Mercury, Copper, Cadmium and Manganese Hollow Cathode lamps are supplied with the system.
Flame analysis system:	
Sensitivity	Cu 5µ/ml > 0.850 Abs.
Burner Head	Titanium alloy burner
Nebuliser	High efficiency glass nebuliser
Atomisation chamber	Corrosion resistant material
Position adjustment	Automatic flame burner optimum height setting via computer.
Safety	Automatic ignit on of mixing air-acetylene gas with protection control via computer
Background correction system:	
Deuterium Lamp Background Correction	1.0 Abs
Self reversal Background Correction	1.0 Abs
Data processing:	
Analytical methods available	Flame AA, graphite furnace, hydride and flame emission.
Determination Method	Calibration curves (1-3 order) standard addition method.
Repetition Times	1-20 times with calculations if average value, standard deviation and relative standard deviation.
Result Printout	Output of parameters, data, spectra and calibration curves.
Power supply:	
Power supply	220v 50Hz AC 200W
Hollow cathode Lamp system	2 lamp simultaneously lit (one with pre-heat) Mercury Hollow Cathode Lamp is supplied with the system
Capacitive Discharge Electro-thermal Atomiser complete with furnace head computer power supply, interface and software.	
The integration of flame burner and graphite furnace in one compact unit enables automatic interchanging between the flame mode and furnace mode using few simple keystrokes.	
Uniformity of the graphite tube temperature reduces chemical interferences and residue interference.	
The PC controlled heating program offers a multitude of heating method choices. The optical temperature control ensures rapid heating and extends the life of the graphite tube.	
The furnace will only be activated for heating when the cooling water flow is sufficient to cool down the furnace body.	
Graphite Tube Fracture Protection: Power supply used for heating is automatically shut off and an alarm is given if a graphite tube fracture is detected.	
Characterictic mass	0.8Pg
Temperature Range	Ambient - 2750°C
Heating	Voltage feedback control when drying and ashing. Optical temperature control when atomizing.
Heating Program	Up to 9 steps with choice of ramp, ladder and full power heating.



Atomic Absorption Spectrophotometer

Fully Automatic Combined Flame and Graphite Furnace System

Model AAS-3900

Standard Configuration	
QTY	DESCRIPTION
1	Communication Lead
1	AA/Win Software
2	Nebulisers (one spare)
10	Pyrolytically Coated Graphite Tubes with Platform
1	HCL Hg for wavelength calibration
3	HCL Mn x1, Cu x1, Cd x1
8	HCL Holders
1	Power Lead (US standard)
1	Power Lead (EU standard)
1	Drain Trap and Tubing
1	Pipette & Tips
N/A	Gas Tubing (Various)
N/A	Fittings (Various)
N/A	Tools



Atomic Absorption Spectrophotometer

Fully Automatic Combined Flame and Graphite Furnace System

Model AAS-3900

ADVANCED GRAPHITE FURNACE:

The unique design of the graphite furnace reduces the chemical interference effects and memory effects by uniformly heating the graphite electrode.

The computer controlled heating program allows the user to select the best heating program for the analysis.

The optical temperature during the atomization stage ensures the rapid heating and rapid analysis. This helps to extend the life of the graphite tube and enhances analytical accuracy.

PROVEN SAFETY FEATURES:


The flame conditions are continuously monitored and should the flow rates change an audible alarm sounds.

The pressure of the support gas (oxidant) is monitored constantly. If the pressure changes then the flow of the fuel gas will be stopped and the flame will be safely extinguished.

A sensor monitors the level of liquid in the drain and will prevent ignition if too low. The flame will also be extinguished if the level of liquid in the drain changes significantly.

The argon pressure for the graphite furnace is constantly monitored and should it change the heating cycle for the graphite electrode will immediately cease and the graphite electrode will be de-energized.

Cooling water flow rates for the graphite furnace are also monitored for changes and should changes occur the heating program will cease. If the graphite tube should fracture during the heating program the heating will cease.





Atomic Absorption Spectrophotometer

Fully Automatic Combined Flame and Graphite Furnace System

Model AAS-3900

Accessories

HYDRIDE AND COLD VAPOUR ANALYSIS

- A hydride generator is available for the determination of elements such as Arsenic, Selenium, Antimony, Tellurium and mercury at ultra low levels. The hydride generator is supplied with an absorption cell, and electrical absorption cell heater and controller and all necessary burner fittings.



AAS FLAME AUTOSAMPLER

- Sequential auto-sampler allows the automated analysis of 50 or more samples and calibration standards. The system allows for automatic update of standard values and curve parameters by using up to 8 standards, blanks and QC standards. A double wash station with facility for use of sample blank or pure water for probe wash avoids sample and standard contamination. An inert Teflon probe is supplied.



GRAPHITE FURNACE AUTOSAMPLER

- The graphite furnace auto-sampler system allows automatic update of calibration data, matrix modifications and automatic dilution of samples.






Atomic Absorption Spectrophotometer

Fully Automatic Combined Flame and Graphite Furnace System

Model AAS-3900

Software Specifications

AA-Win Pro Software is a powerful and intuitive software product designed to allow control and data acquisition from the SpectroAAS-3900 Atomic Absorption Spectrophotometer. The AA-Win Pro software allows the Analyst to control all aspects of their analytical method whilst providing an extensive range of tools for data collection, storage and interpretation.




Lamp turret setup, operating and warm-up currents, along with the desired analytical wavelengths are easily selected in this configuration.



No.	Wavelength	Sample	Wavelength	Wavelength	Wavelength	Wavelength	Wavelength	Wavelength	Wavelength
1									
2									
3									
4									
5									
6									
7									
8									
9									
10									
11									
12									
13									
14									
15									
16									
17									
18									
19									
20									
21									
22									
23									
24									
25									
26									
27									
28									
29									
30									
31									
32									
33									
34									
35									
36									
37									
38									
39									
40									
41									
42									
43									
44									
45									
46									
47									
48									
49									
50									

Use the sample table to perform quick measurements of both Standards and Samples. Easily append the sample table to add new samples or even revise calibration curves either by manual introduction or using an Autosampler.

Ensure optimal peak position at the chosen analytical line by scanning the emissions spectra.




Each stage of analysis setup is made quick and simple by means of the Sample Wizard.



View up-to-date calibration curves in 1st, 2nd, or 3rd order using a standard calibration or standard addition. Perform retrospective curve fits to ensure optimum correlation.



Obtain reliable and accurate results by using the Energy control feature to manually optimize atomiser position and setup, Use the Auto-balance feature to ensure energy level, and optical alignments are optimized when using background correction.



View real-time signal acquisition for flame, graphite furnace and hydride generation analysis.



Atomic Absorption Spectrophotometer

Fully Automatic Combined Flame and Graphite Furnace System

Model AAS-3900

Periodic table of chemical elements | Details

1A	2A	3B	4B	5B	6B	7B	8	9	10	11	12	13	14	15	16	17	18
H												He					
Li	Be											B	C	N	O	F	Ne
Na	Mg											Al	Si	P	S	Cl	Ar
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	Kr
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	Xe
Cs	Ba	La	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At	Rn
Fr	Ra	Ac															

Legend:
Flame: (white box)
Furnace: (blue box)
Hydride: (purple box)

Legend elements:
Flame: Ce, Pr, Nd, Pm, Sm, Eu, Gd, Tb, Dy, Ho, Er, Tm, Yb, Lu
Furnace: Th, Pa, U, Np, Pu, Am, Cm, Bk, Cf, Es, Fm, Md, No, Lr
Hydride: (none listed)

Elements that can be analyzed with the AAS-3900

Software

- The user friendly software requires a Windows platform and operates within Win95, Win98, Win NT, Win 2000 and WinXP. The system uses a number of software wizards to guide the operator through setting up procedures.
- The software controls the automatic switch over for the Hollow Cathode Lamps and automatically optimizes working parameters for the system. The software also allows manual input of data to ensure that the operator always stays in control. The software will automatically complete the configuration of the system for analysis.
- The user has the choice of two methods of background correction namely the self reversal system OR the traditional deuterium lamp background correction system.
- During the analysis cycle of both the flame and graphite furnace the software shows the entire measurement process. This includes measured values, temperature steps, time etc. all signal and temperature data is stored for future re-call and printout.
- Detailed reporting and QC control software is included within AAWin allows printout of spectra, standard calibration curves, analysis and signal data. Full printout of operating parameters is also available for user references.
- The following methods of analysis can be carried out using the AAS-3900 system Absorption, emission, graphite furnace analysis, hydride and cold vapour analysis.

