


Grip Strength Meter for Mice and Rats

To a ssess m uscle hyperalgesia this i nstrument measures the forelimb grip force u sing a digital force transducer.

The transducer is connected to a wire mesh grid connected to an anodized base plate. The animal is held by its tail and is gently passed over the mesh until it grasps the grid with its front paws. Three grip force measures are made and the average of these readings represents the animal's forelimb grip force at that particular time.

The last maximum force is held in a "peak & hold" type readout until reset. The unit gives readings in one gram increments. Both rat and mouse grids are supplied and can be easily exchanged. Maximum force range is 2000 Gr (10% over range allowed). Higher ranges available on special order.

The heavy base plate has suction feet to resist even large pulling forces. Both battery and line adapter operation is possible.





Part #	Description
2200	Grip Strength Meter for Mice and Rats
2250	Grip Strength Stand and Meshs Only