

Technical specification

Cherish your life, Cherish your health!

Ventilation mode

IPPV, V-A/C, P-A/C, P-SIMV, V-SIMV, PCV, PSV, SPONT/CPAP, PRVC(optional), APRV(optional), Duol level(optional), SIGH, MANU

Ventilator parameter

Tidal volume (Vt) 0, 20 ~ 2000 ml

Frequency (Freq) 1 min ~ 100 min

Oxygen concentration 21 % ~ 100 %

I:E 4: 1 ~ 1: 8

PEEP 0 cmH₂O ~ 40 cmH₂O

Pressure limit 20 cmH₂O ~ 100 cmH₂O

Pressure triggering sensitivity (Ptr) -20 cmH₂O ~ 0 cmH₂O (Based on PEEP)

Flow triggering sensitivity (Ftr) 0.5 cmH₂O ~ 30 cmH₂O

SIGH On/off 1/100 ~ 5/100

Apnea ventilation OFF, 5 s ~ 60 s

Monitoring parameter

Frequency (Freq) 0/min ~ 100/min

Tidal volume (Vt) 0 ml ~ 2500 ml

MV 0 L/min ~ 99 L/min

Airway pressure 0 cmH₂O ~ 100 cmH₂O

Dynamic lung compliance monitoring 1 mL/cmH₂O ~ 1000 mL/cmH₂O

Oxygen concentration 15 % ~ 100 %

Packing size

Main engine: L 560 * W 560 * H 605 mm

G.W.: 35 KG, N.W.: 17 KG

Air compressor: L 683 * W 687 * H 1140 mm

G.W.: 100 KG, N.W.: 65 KG

Alarm and protection

The AC power failure alarm Power failure or no connection

Internal backup battery low voltage alarm 11.3 ± 0.3 V

No tidal volume No tidal volume within 6 s

High minute volume alarm 5 L/min ~ 99 L/min

Low minute volume alarm 1 L/min ~ 30 L/min

High airway pressure alarm 20 cmH₂O ~ 100 cmH₂O

Low airway pressure alarm 0 cmH₂O ~ 20 cmH₂O

High oxygen concentration alarm 19 % ~ 100 %

Low oxygen concentration alarm 18 % ~ 99 %

Continuous pressure alarm (PEEP + 1.5 cmH₂O) over 16s

Suffocation warning 5 ~ 60 s

Fan error Show on screen

Oxygen deficit Show on screen

The maximum limited pressure < 12.5 kPa

Working condition

Gas source O₂ Air

Pressure 280 kPa - 600 kPa

Voltage 220 V ± 22 V

Power frequency 50 Hz ± 1 Hz

Input power 900 VA (With air compressor)

250 VA (Without air compressor)

Oscillogram

P-T (Pressure-Time)

F-T (Flow-Time)

P-V Loop (Pressure-Volume Loop)

S1100 ICU Ventilator

Friendly Powerful Reliable





S1100 ICU Ventilator

Application

The ventilator makes a good performance in operation room, ICU department and emergency treatment. It is used to assist or replace the spontaneous breathing for adult, pediatric and neonatal more than 2kg. 25 years experience in market-oriented ventilator make us professional and reliable, satisfying all your needs in ventilation. Due to the flexible configuration, good quality and competitive price, S1100 has soon become the superstar of the market.


Feature

- 15" TFT screen displays the ventilation parameters, alarm information, and oscillograms, make every operation more easily
- Multiple ventilation mode can meet different clinical requirements. (IPPV A/C PCV SIMV PSV SPONT/CPAP SIGH MANUAL)
- 3 oscillograms for your choice, 2 of them can be displayed on the screen at the same time
- Humidifier can heat and wet breathing gas, makes it comfortable for patient to breathe
- Rapid oxygen supply, automatically offer high flow rate oxygen within two minutes
- High temperature resistance breathing circuit is reusable and anti-pollution

Safety

- 14 types of sound and visual alarm information, easier for users to do some error checking and troubleshooting
- Built-in oxygen concentration sensor, ensure stable precision of oxygen concentration
- Easy to move with four casters, easy to stop with two brakes
- Separate design of electronic circuit and gas flow rate keep safe running of ventilator
- Compact long life internal battery can provide emergency power, avoid risk of patient
- Self-check before operation, eliminate system mistake

Friendly Powerful Reliable



Packing list

Name: Ventilator

Type: S1100

NO.	Name	Specification	Quantity	Unit	Remark
1.	Main body	2.782.075A	1	PC.	Packing box 1
2.	Power cord	1.2M	1	PC.	
3.	Air compressor	2.782.134A	1	PC.	Packing box 2
4.	High pressure hose for air compressor	4.472.023	1	PC.	
5.	High pressure O ₂ hose	4.472.001	1	PC.	
6.	CO ₂ sensor	KM7012	1	PC.	
7.	Power cord	4M	1	PC.	
8.	Fuse	T2A	4	PCS.	
9.	Fuse	T4A	2	PCS.	
10.	Fuse	T10A	2	PCS.	
11.	Fuse	T6.3A	2	PCS.	
12.	Screw	M5X16	6	PCS.	
13.	Serrated saddle	φ 5	6	PCS.	
14.	Hanger		1	PC.	
15.	Humidifier (adult)	SH330	1	PC.	
16.	Breathing tube (C- I -1)	(φ 22- φ 22) × 800mm	2	PCS.	
17.	Breathing tube(C- I -1)	(φ 22- φ 22) × 400mm	3	PCS.	
18.	Breathing tube(C- I -2)	(φ 22- φ 22) × 400mm	3	PCS.	
19.	Breathing tube(C- I -2)	(φ 22- φ 15) × 800mm	2	PCS.	
20.	Water trap	Adult	2	PCS.	
21.	Y connector	Adult	1	PC.	
22.	Y connector	Child	1	PC.	
23.	L connector		2	PCS.	
24.	Test lung	1L	1	PC.	
25.	Operation manual		1	PC.	
26.	Certificate of Quality		1	PC.	
27.	Packing list		1	PC.	

Main technical specification of S1100

Power supply

- Voltage 100-240 V ~
- Frequency 50 Hz/60 Hz
- Input power 80 VA
- Internal power DC12V 4Ah
- Fuse
 - Main power T2AH 250V
 - Internal power T4AH 250V

Ventilation mode

- IPPV
- A/C
- V-SIMV
- P-SIMV
- PSV
- PCV
- SPONT/CPAP
- PRVC
- APRV
- DuoLevel
- SIGH
- MANUAL

Ventilation parameter adjustment

- Frequency
 - Adjusting range 1 /min ~ 100 /min
(Under SIMV: 1 /min ~ 40 /min
All mode except SIMV: VT_H 4 /min ~ 40 /min
 VT_L 20/min ~ 100 /min)
 - Allowable error ± 1 /min or \pm (10% setting value), whichever is the greater
- Inhalation time (I:E) (T_{insp})
 - Adjusting range 0s ~ 12s(except SIMV mode, I:E 4:1 ~ 1:8)
 - Allowable error ± 15 %
- Tidal volume (V_T)
 - Adjusting range 0 mL ~ 1500 mL
 - And:
 - VT_H 250 mL ~ 1500 mL
 - VT_L 0. 20 mL ~ 300 mL

- Allowable error $\pm (10 \text{ mL} + 10\% \text{ setting value})$
- Minute ventilation (MV)
 - Max MV
 - $VT_H > 18 \text{ L/min}$
 - $VT_L > 10 \text{ L/min}$
- PEEP
 - Adjusting range $0 \text{ cmH}_2\text{O} \sim 20 \text{ cmH}_2\text{O}$
 - Allowable error $\pm (2 \text{ cmH}_2\text{O} + 5\% \text{ setting value})$
- CPAP
 - Adjusting range $0 \text{ cmH}_2\text{O} \sim 20 \text{ cmH}_2\text{O}$
 - Allowable error $\pm (2 \text{ cmH}_2\text{O} + 5\% \text{ setting value})$
- Continuous flow
 - Adjusting range $VT_H 7 \text{ L/min} \sim 60 \text{ L/min}$
 $VT_L 2 \text{ L/min} \sim 30 \text{ L/min}$
 - Allowable error $> 3\text{L/min}, \pm 15\%$
 $< 3\text{L/min}, \pm 0.5 \text{ L/min}$
- Inspiratory trigger
 - Pressure trigger (PTr)
 - Adjusting range $- 20 \text{ cmH}_2\text{O} \sim 0 \text{ cmH}_2\text{O} \text{ (Under PEEP)}$
 - Allowable error $\pm (1 \text{ cmH}_2\text{O} + 10 \% \text{ setting value)}$
 - Flow trigger
 - Adjusting range $0.5\text{L/min} \sim 30 \text{ L/min}$
 - Allowable error $> 3\text{L/min}, \pm 15\%$
 $< 3\text{L/min}, \pm 0.5 \text{ L/min}$
- Pressure control (Pc)
 - Adjusting range $5 \text{ cmH}_2\text{O} \sim 60 \text{ cmH}_2\text{O}$
 - Allowable error $\pm (2 \text{ cmH}_2\text{O} + 5\% \text{ setting value})$
- Pressure support (Ps)
 - Adjusting range $0 \text{ cmH}_2\text{O} \sim 60 \text{ cmH}_2\text{O}$
 - Allowable error $\pm (2 \text{ cmH}_2\text{O} + 5\% \text{ setting value})$
- O₂ concentration
 - Adjusting range $21 \% \sim 100 \%$
 - Allowable error $\pm 3 \% \text{ (v/v) or } \pm 10 \% \text{ setting value, whichever is the greater}$
 - Response time of the VENTILATOR from a volume fraction of 21 % to 90 %
 - $V_t 500 \text{ mL} \quad 2\text{s}$
 - $V_t 150 \text{ mL} \quad 2.5 \text{ s}$
 - $V_t 30 \text{ mL} \quad 3\text{s}$
- Holding time (Inhalation platform) (TIP)
 - Adjusting range $0 \sim 6\text{s}(0\% \sim 50\% \text{ inhalation time})$

allowable error $\pm 0.1s$ or $\pm (10\% \text{ set value})$, whichever is the greater

- Apnea ventilation
 - Adjusting range OFF. 5s - 60s
- Maximal Inspiratory Flow Rate
 - Maximal Inspiratory Flow Rate: $> 60 \text{ Umin}$
- Max set working pressure (pressure limit range) 20cmHP - 100cmH₂O
- Max pressure $< 125\text{cmH}_2\text{O}$
- Output MV for manual ventilation $> 2\text{SL}$ min
- Nebulizer gas maximum output pressure $< 0.2 \text{ MPa}$
maximum output flow $> 8 \text{ Umin}$

explanation:

- 1) set pressure limit value to produce maximum working pressure (the upper limit of high airway pressure alarm)
- 2) do not use positive pressure at exhalation
- 3) maximum limit pressure is pressure of safety valve.

Ventilation parameter Monitoring range

- The following are monitoring parameters under the environment of body temperature and pressure saturated
- Fig display:


Parameter	Range	Resolving power	Accuracy
Frequency (Freq)	0/min - 100/min	1/min	$\pm 1/\text{min}$ or $\pm (10\% \text{ actual range})$, whichever is the greater
Tidal volume (VT)	0 mL - 2000mL	10ml	$> 50\text{mL}$, $\pm (4.0 \text{ mL} + 1.5\% \text{ actual range})$ ($< 50 \text{ ml}$, $\pm 1 \text{ mL}$)
Minute ventilation (MV)	0 Umin - 99 Umin	0.1 Umin	$> 3 \text{ Umin}$, $\pm (15\% \text{ actual range})$ ($< 3 \text{ Umin}$, $\pm 0.5 \text{ Umin}$)
Peak airway pressure (P _{peak})	0 cmH ₂ O - 100 cmH ₂ O	1 cmH ₂ O	$\pm (2\% \text{ full scale} + 4\% \text{ actual range})$
O ₂ concentration	15% - 100%	0.1 %	$\pm [2.5\% (v/v) + 2.5\% (O_2 \text{ concentration level})]$
Compliance	1~1000 mL/cmH ₂ O	1 mL/cmH ₂ O	

- Waveform display:
 - Time- airway pressure (under all modes)
 - Time—flow (under all modes)
 - pressure volume loop (all modes)

Remarks: The machine cannot record all the adjusting and monitoring values.

- the purpose, sensor position, type and sampling method of control, measurement and display device

the purpose	sensor position	type	sampling method
Airway pressure	Exhaling end	pressure- voltage、 Simulated data	Choose average value from multiple sampling
PEEP	Exhaling end	pressure- voltage、 Simulated data	Choose average value from multiple sampling
Continuous pressure	Exhaling end	pressure- voltage、 Simulated data	Choose average value from multiple sampling
freq	built-in	time、 simulated data	Choose average value from multiple sampling
Inhalation time	built-in	time、 simulated data	Choose average value from multiple sampling
Holding time	built-in	time、 simulated data	Choose average value from multiple sampling
tidal volume	Exhaling end	flow- voltage、 Simulated data	Choose average value from multiple sampling
flow	Exhaling end、 inhaling end	flow-voltage、 simulated data	Choose average value from multiple sampling
O ₂ concentration	Inhaling end	O ₂ concentration-voltage、 simulated data	Choose average value from multiple sampling



编号 120191 的应 O i:1071 20 J ti9

若 业执 照 闺 本



门批... 礼 讲<N*- 抽
1: 讲,北 讲 井 if @CI, : 力,
并 次 (. 更 名 妒 己
而 . 讲 可 队 肖 佑 息.

统一社会信用代码

913201166673974138 (1/1)

名称 南京舒普思达医疗设备有限公司
类型 有限责任公司

注册资本 1921.2(?.万元整)

成立日期 2008年01月18日

法定代表人 楼田甘甘

营业期限 2008年01月18日至2028年01月17日

经营范围 医疗器械生产（按许可证所列项目生产、经营），从事医疗器械研发；自产和代理各类商品及技术的进出口业务（国家限制企业经营或禁止进出口的商品和技术除外）；仪器仪表的研发、生产、销售；动物用麻醉、呼吸、监护、养护设备的研发、生产、销售；机械电子设备、试验设备制造、加工；如生产、材料；软硬件及工业设计技术咨询、技术服务。（依法须经批准的项目，经相关部门批准后方可开展诊疗活动）

住 所 南京市江化浦区溧江路0号6楼 压 寸县

登记机关 2019 年 07 月 12 日



医疗器械 产、 可证

许可证编号：苏食药监械生产许20070089号

企业名称：

法定代表人：楼[]甜

：产产品 抒记 表

企业负责人：彭家良

住 所：南京中江北新区冈富路9号6号
楼[]层、[]次

发证部门：江 苏

有效期限：至 2022 年 。1 月 。8 日

发证日期： 2017

0 9 日

上 许

南京舒普思达医疗设备有限公司



该复印件与原件一致
仅供存档备案，不做经营凭证
再复印无效
生产范围：
公司名称

生产地址：南京市江北
层、三层

生产范围：见医疗器械

省药

• Ridlerstrane 65·80339 Munich • Germany

