

# YX18(24) LDJ

Portable Stainless Steel Steam Sterilizer



- The product uses saturated steam to rapidly and efficiently sterilize articles.
- It can be used in hospital, public health center, health stations, clinics in factories and mines, etc. to sterilize medical equipment, surgical dressing, glass utensils, solutions, substrata, etc.
- It can also be used as a cooker in plateau areas and be used to make high-quality drinking water in industrial and mining enterprises.
- Auto power cut-off of overheating protection for water-off & Current overload.

## SPECIFICATION

Volume of the disinfecting pail	:	- 0.018 m <sup>3</sup> (Ø 280 x 265 mm) - 0.024 m <sup>3</sup> (Ø 280 x 265mm)
Max. Working Pressure	:	0.165 Mpa
Temperature	:	126°C
Working Pressure	:	0.14 Mpa
Mpa Adjustment	:	105°C - 126°C
Timer	:	0 - 60 min
Heat Average	:	≤ ± 1°C
Power Supply	:	AC 220 V, 50 Hz
Heater Power	:	1.5 KW
Net Weight	:	14 kg