



CE



BT-LED4D

Shadowless Operating Lamp

Zhangjiagang Bestran Technology Co.,Ltd






Precise dimming control for various lighting request

With 16-stage illuminance adjustment, 4-stage color temperature adjustment, 16-stage spot size adjustment by electrical control, endo mode, R9 mode, green light (GL) mode.

Unique maintenance room for quick and easy maintenance

The innovative design of photoelectric separation keeps the circuit board away from heat sources. It increases the service life time, also improves the simplicity, economy and timeliness of the maintenance. What's more, the control box can also be easily opened and maintained.





BESTRAN

BT-LED4D

Shadowless Operating Lamp

Model: BT-LED4D

Color: White + dark gray

LED brand:OSRAM

Bulb qty.(pcs):112PCS

Color temperature:

1st-stage 3500K 2nd-stage 4000K

3rd-stage 4500K 4th-stage 5000K

Illumination:160,000 Lx

CRI(Ra): ≥ 95 , Peak 96; R9: 92; R13: 95

Illuminance depth:1000mm