



UV PASS BOX

WITHOUT HEPA FILTERS

LUP-B13, LUP-B14, LUP-B15

UV Pass Box (without HEPA filters) LUP-B13, LUP-B14, LUP-B15 are static type UV pass box equipped UV germicidal lamp, both side electrical interlocked doors, and high quality see through glass window. Pass box is made up of durable high quality 304 grade stainless steel, and updated quality parts to meet the demands of controlled clean room environment operations.

Features

- Indicator lights to know the door lock/unlock condition
- Electronically interlocked door type
- Red press button for on-off operations
- Certified with CE and ISO standards

Applications

Used to transfer articles, materials between different levels of clean room across pharmaceutical, food, life science and other research, quality inspection industries.

Specification

Model no.	LUP-B13	LUP-B14	LUP-B15
External Dimension (W*D*H)	940 × 860 × 880 mm	1040 × 960 × 980 mm	1140 × 1060 × 1080 mm
Internal Dimension (W*D*H)	800 × 800 × 800 mm	900 × 900 × 900 mm	1000 × 1000 × 1000 mm
UV Lamp	8 to 11 W × 1 lamp		
Door Interlock Type	Electronical		
Material	304 Stainless Steel		
Net Weight	45 kg	50 kg	55 kg
Packaging Dimension (W*D*H)	990 × 900 × 980 mm	1090 × 1010 × 1080 mm	1190 × 1110 × 1180 mm
Power Supply	220 V, 50 Hz		