

Haier

-86°C ULT Freezer
DW-86W100J



DW-86W100J



Qingdao Haier Biomedical Co., Ltd.

Model		DW-86W100J	
Technical Data	Cabinet Type	Chest	
	Climate Class	N	
	Cooling Type	Direct cooling	
	Defrost Mode	Manual	
	Refrigerant	HC	
	Sound level(dB(A))	49	
Performance	Cooling Performance(°C)	-86	
	Temperature Range(°C)	-40~-86	
Control	Controller	Microprocessor	
	Display	LED	
Electrical Data	Power Supply(V/HZ)	220~240/50	
	Power(W)	550	
	Electrical Current(A)	4	
Dimensions	Capacity(L/Cu.Ft)	100/3.5	
	Net/Gross Weight(approx)	kg	138/160
		lbs	304.2/352.7
	Interior Dimension(W*D*H)	mm	470*450*480
		in	18.5*17.7*18.9
	Exterior Dimension(W*D*H)	mm	769*825*1120
		in	30.3*32.5*44.1
	Packing Dimension(W*D*H)	mm	845*855*1250
in		33.3*33.7*49.2	
Container load(20'/40'/40'H)	12/24/48		
Accessories	Caster	Y	
	Foot	Y	
	Porthole	Y/1	
	Shelves/Inner doors	-/1	
	USB Interface	Y	
	Remote Alarm	Y	
	5V Power Supply Port	N/A	
	Temperature Recorder	Optional	
	Rs232/485 Port	Optional	
	CO ₂ Backup System	Optional	
	LN ₂ Backup System	Optional	
Other	Certification	CE	

Distributed by:



Qingdao Haier Biomedical Co., Ltd.

Haier industrial Park No 1 Haier Road, Qingdao,
266101, China