



Technical Data & Ordering Information

Model Name	CKD015	CKD065
No. of positions	1	6
Flask Capacity	500 ml Kjeldahl Flask. Also compatible with 800 ml Kjeldahl Flask	
Element Temperature	450°C	
Temperature regulators	1	6
Display	LED Display with time, safe temperature & process status	
Timer	1-999 Minutes	
Power consumption	150W	900W
Electrical regulation	0.65A	3.9A
HSN Code:	84198990	
Ordering Code	BLFACKD01500000000	BLFACKD06500000000

Accessories

Model Name	Ordering Code
Graham Condenser 200mm Jacket ACCS	BLA0BGW200KJDLGHM0
CLASSICAL KJELDHAL SPLASH HEAD-2429 ACCS	BLA0KSH24291926CKJ
CKD065 without Stand ASM	BLAACKD06500000001
CKD065 with Stand ASM	BLAACKD06500000000
CKD015 Stand ASM	BLAHCKD01510000000

Note : All accessories are provided by default with the unit