

Double door fully automatic Washer disinfector

Application: surgical instruments

Automatic Glassware Washer



WD-220M/-320M



WD-420M

Advantage

1. Front and back two toughened glass doors, prevent cross infection.
2. Efficiently removes organic residues and metal ions through automated high-temperature, high-pressure water cleaning.
3. Compact and lightweight design with corrosion-resistant stainless steel construction, suitable for small laboratories.
4. Microcomputer control system with LCD screen for easy operation and customization of cleaning modes and parameters.
5. Efficient cleaning with a powerful imported pump, even water distribution, secure sealing, and rotating spray arms for thorough cleaning.
6. Secure door lock with emergency opening, comprehensive protection with alarms, and English prompts for safe operation.

Specification

Model	WD-220M	WD-320M	WD-420M
Capacity(L)	220	320	420
Voltage(V/Hz)	380V/50~60HZ;220V~240V/50~60HZ		
Washer Load Level	1	1	1
Water Wash Power(Kw)	0.75	1.5	1.5
Water Heater Power(Kw)	15		
Dry Power(Kw)	3.6		
MAX. Power(Kw)	17		
Water Consumption/Cycle(L)	22	24	26
Tank internal temperature (C)	RT+5~99		
Noise(dBa)	<60		
Water Supply Pressure (Mpa)	0.3-0.8		
Chamber Size(WxDxH mm)	600x629x658	600x629x820	688x700x703
External Size(WxDxH mm)	960x800x1880	960x800x2030	1100x900x2165
Shipping Size (WxDxH mm)	1200x1000x2200	1200x1000x2350	1300x1040x2350
N.W./G.W. (kg)	270/360	290/370	370/460

Cleaning Rack



Specification

Cleaning Rack for Glassware Washer (except for WD-60 and WD-200I)

UniversaRacks			
Model	CWIR30	CWIR36	CWIR60
Picture			
Description	To wash Max. 30ea of 250ml infusion bottles and caps; or other volumetric flasks conical asks or beakers	To wash Max.36ea of 10~250ml volumetric flasks, conical flasks or beekers(diameter within 90mm)	To wash Max.60ea of 10~250ml volumetric flasks, conical flasks or beekers(diameter within 60mm)
Model	CWIR120	CWIR222	CWIR260
Picture			
Description	To wash 120 ea of centrifuge tubes, test tubes, 50ml volumetrics flask	222position of injection tubes and adjustable trolleies. It can wash 222ea of tubes or vials and capes.	260 washable injection vials

Cleaning Rack

Surgical Instruments

Model		
Picture		
Description	K-05 for surgical scissors	K-04 1/4 Fan-shaped Test Tube Cleaning Basket Bottle volume: 10 - 20ml Clean test tubes, centrifuge tubes, etc. Four 1/4 modules form a layer of cleaning basket
		K-02 Large Quantity Cleaning Frame Mmodule Bottle volume: 10 - 2000ml Clean jars, glass rods, bottle caps Two 1/2 modules form a layer of cleaning basket

Petri Dish Cleaning Rack

Model	CWPDR9070L	CWPDR9070U
Picture		
Description	70 petri dishes with diameter less than 90mm (height less than 30mm) can be cleaned. This model is limited to the upper layer.	75 petri dishes with diameter less than 100mm (height less than 30mm) can be cleaned. This model is limited to lower layer.