

Gel Card Centrifuge

Application:

Gel card centrifuge is used in blood serology, blood type routine test, erythrocyte washing, microcolumn agglutination, immunology and other tests.



Features:

- * DC Brushless motor, low noise.
- * Microprocessor control system & soft touch control panel.
- * LCD display: RCF, Speed, Time.
- * Electric lid interlock, over speed and imbalance protection.
- * 20 programs for standard and programmability are available.

Technical Parameters:

| Model | BKC-GC12 | BKC-GC24 |
|------------------------|-------------------------|-------------|
| Max. Speed | 4000rpm | |
| Max. RCF | 2250xg | |
| Rotor | 12 cards | 24 cards |
| Speed Precision | ±20rpm | |
| Time Range | 0~99min | |
| Noise | ≤55dB | |
| Power Supply | AC110/220V±10%, 50/60Hz | |
| External Size(W*D*H)mm | 460*320*230 | 500*360*250 |
| Package Size(W*D*H)mm | 560*395*330 | 610*435*340 |
| Gross Weight(kg) | 18 | 27 |

Table Top Dairy Centrifuge



Application:

It is special equipment for milk and dairy analysis, such as analyze the samples after fat separation through pasteurization, Gerber's method, etc.

Features:

- * The LCD display indicates the speed, time and temperature.
- * Microprocessor control.
- * Electric lid lock, over speed, over temperature and imbalance protection.
- * Stainless steel centrifugal chamber and lid.

Technical Parameters:

| Model | BKC-MF5A |
|------------------------|---------------------|
| Max. Capacity | 8*Gerber tubes |
| Max. Speed | 1500rpm |
| Max. RCF | 420xg |
| Speed Precision | ±50rpm |
| Temp. Range | RT+10°C~65°C |
| Timing Range | 0~99h59min |
| Noise | ≤60dB |
| Consumption | 1.5KW |
| Power Supply | AC110/220V, 50/60Hz |
| External Size(W*D*H)mm | 520*415*330 |
| Package Size(W*D*H)mm | 720*610*530 |
| Net/Gross Weight(kg) | 35/60 |

Rotors for BKC-MF5A:

| Rotor No. | Rotor | Capacity | Speed(rpm) | RCF(×g) |
|-----------|-------------|----------------|------------|---------|
| No.1 | Angle Rotor | 8*Gerber tubes | 1500 | 420 |